

Application

- Vascular Surgery
- Flushing of veins

Technical specification

- Pack of 3 - two hybrids (one large and one small) and one 6fg rounded end fine Cannula
- Hybrid version of the traditional cone and olive – cone to accommodate a wide range of sizes and olive with curved ‘waist’ for easy grip on the vein
- Standard luer connection
- Designed to match quality of reusable equivalent
- 25 sets per box

Additional Features and Benefits

- Non blocking design
- Clean room moulded to eliminate manufacturing debris
- Supplied sterile in an individual peel pouch

